

NARF Antibody



PACO03845

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

NARF

Uniprot

Q9UHQ1

Synonyms:

NARF; Nuclear prelamina A recognition factor; Iron-only hydrogenase-like protein 2; IOP2

Immunogen:

Synthesized peptide derived from the Internal region of human NARF.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A