

Phospho-GFAP (S38) Antibody



PACO03855

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

GFAP

Uniprot

P14136

Synonyms:

GFAP; Glial fibrillary acidic protein; GFAP

Immunogen:

Synthesized peptide derived from human GFAP around the phosphorylation site of S38.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A