

GJA3 Antibody



PACO03869

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

GJA3

Uniprot

Q9Y6H8

Synonyms:

GJA3; Gap junction alpha-3 protein; Connexin-46; Cx46

Immunogen:

Synthesized peptide derived from the Internal region of human Connexin 46.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A