

AMPD3 Antibody



PACO03885

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

AMPD3

Uniprot

Q01432

Synonyms:

AMPD3; AMP deaminase 3; AMP deaminase isoform E; Erythrocyte AMP deaminase

Immunogen:

Synthesized peptide derived from the Internal region of human AMPD3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A