

ARFIP1 Antibody



PACO03886

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

ARFIP1

Uniprot

P53367

Synonyms:

ARFIP1; Arfaptin-1; ADP-ribosylation factor-interacting protein 1

Immunogen:

Synthesized peptide derived from the C-terminal region of human Arfaptin 1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A