

ARHGEF16 Antibody



PACO03887

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

ARHGEF16

Uniprot

Q5VV41

Synonyms:

ARHGEF16; EPHEXIN4; NBR; Rho guanine nucleotide exchange factor 16; Ephexin-4

Immunogen:

Synthesized peptide derived from the Internal region of human ARHGEF16.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A