

# GNRH1 Antibody



PACO03909

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

IHC:1:100-1:300

**Protein Background:****Gene ID:**

GNRH1

**Uniprot**

P01148

**Synonyms:**

GNRH1; GNRH; GRH; LHRH; Progonadoliberin-1; Progonadoliberin I

**Immunogen:**

Synthesized peptide derived from the Internal region of human GnRH I.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A