

OR4K3 Antibody



PACO03923

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR4K3

Uniprot

Q96R72

Synonyms:

OR4K3; OR4K3P; Olfactory receptor 4K3; Olfactory receptor OR14-14

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4K3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A