

# EPHA10 Antibody



PACO03938

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

EPHA10

**Uniprot**

Q5JZY3

**Synonyms:**

EPHA10; Ephrin type-A receptor 10

**Immunogen:**

Synthesized peptide derived from the Internal region of human EphA10.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A