

GPR25 Antibody



PACO03939

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

GPR25

Uniprot

O00155

Synonyms:

GPR25; Probable G-protein coupled receptor 25

Immunogen:

Synthesized peptide derived from the Internal region of human GPR25.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A