

# EPHA6 Antibody



PACO03942

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

EPHA6

**Uniprot**

Q9UF33

**Synonyms:**

EPHA6; EHK2; HEK12; Ephrin type-A receptor 6; EPH homology kinase 2; EHK-2; EPH-like kinase 12; EK12

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human EphA6.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A