

DOK7 Antibody



PACO03946

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

DOK7

Uniprot

Q18PE1

Synonyms:

DOK7; C4orf25; Protein Dok-7; Downstream of tyrosine kinase 7

Immunogen:

Synthesized peptide derived from the N-terminal region of human Dok-7.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A