

GPR150 Antibody



PACO03948

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

GPR150

Uniprot

Q8NGU9

Synonyms:

GPR150; Probable G-protein coupled receptor 150

Immunogen:

Synthesized peptide derived from the Internal region of human GPR150.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A