

# GPR42 Antibody



PACO03955

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

GPR42

**Uniprot**

O15529

**Synonyms:**

GPR42; GPR42P; G-protein coupled receptor 42

**Immunogen:**

Synthesized peptide derived from the N-terminal region of human GPR42.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A