

# GRIN2A/GRIN2B Antibody



PACO03978

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

GRIN2A/GRIN2B

**Uniprot**

Q12879/Q13224

**Synonyms:**

GRIN2A; NMDAR2A; Glutamate [NMDA] receptor subunit epsilon-1; N-methyl D-aspartate receptor subtype 2A; NMDAR2A; NR2A; hNR2A; GRIN2B; NMDAR2B; Glutamate [NMDA] receptor subunit epsilon-2; N-methyl D-aspartate receptor subtype 2B; NMDAR2B; N

**Immunogen:**

Synthesized peptide derived from human NMDAepsilon1/2 around the non-phosphorylation site of Y1246/1252.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A