

# GRM8 Antibody



PACO03998

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

GRM8

**Uniprot**

O00222

**Synonyms:**

GRM8; GPRC1H; MGLUR8; Metabotropic glutamate receptor 8; mGluR8

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human mGluR-8.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A