

H2AFZ Antibody



PACO04025

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

H2AFZ

Uniprot

P0C0S5

Synonyms:

H2AFZ; H2AZ; Histone H2A. Z; H2A/z

Immunogen:

Synthesized peptide derived from human Histone H2A around the non-acetylation site of K5.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A