

Phospho-H3F3A (S31) Antibody



PACO04028

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

H3F3A

Uniprot

P84243

Synonyms:

H3F3A; H3.3A; H3F3; PP781; H3F3B; H3.3B; Histone H3.3

Immunogen:

Synthesized peptide derived from human Histone H3.3 around the phosphorylation site of S31.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A