

Phospho-HTT (S421) Antibody



PACO04038

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

HTT

Uniprot

P42858

Synonyms:

HTT; HD; IT15; Huntingtin; Huntington disease protein; HD protein

Immunogen:

Synthesized peptide derived from human Huntingtin around the phosphorylation site of S421.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A