

Phospho-HDAC1 (S421) Antibody



PACO04040

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

HDAC1

Uniprot

Q13547

Synonyms:

HDAC1; RPD3L1; Histone deacetylase 1; HD1

Immunogen:

Synthesized peptide derived from human HDAC1 around the phosphorylation site of S421.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A