

HLA-DPA1 Antibody



PACO04047

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:**Gene ID:**

HLA-DPA1

Uniprot

P20036

Synonyms:

HLA-DPA1; HLA-DP1A; HLABB; HLA class II histocompatibility antigen; DP alpha 1 chain; DP(W3); DP(W4); HLA-SB alpha chain; MHC class II DP3-alpha; MHC class II DPA1

Immunogen:

Synthesized peptide derived from the Internal region of human HLA-DPalpha1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A