

Phospho-APC (S2054) Antibody



PACO04067

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

APC

Uniprot

P25054

Synonyms:

APC; DP2.5; Adenomatous polyposis coli protein; Protein APC; Deleted in polyposis 2.5

Immunogen:

Synthesized peptide derived from human APC around the phosphorylation site of S2054.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A