

AGBL3 Antibody



PACO04099

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

AGBL3

Uniprot

Q8NEM8

Synonyms:

AGBL3; CCP3; Cytosolic carboxypeptidase 3; ATP/GTP-binding protein-like 3

Immunogen:

Synthesized peptide derived from the Internal region of human CCP3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A