

CYP26C1 Antibody



PACO04101

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

CYP26C1

Uniprot

Q6V0L0

Synonyms:

CYP26C1; Cytochrome P450 26C1

Immunogen:

Synthesized peptide derived from the Internal region of human CYP26C1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A