

# OR2T3/OR2T34 Antibody



PACO04112

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR2T3/OR2T34

**Uniprot**

Q8NH03/Q8NGX1

**Synonyms:**

OR2T3; Olfactory receptor 2T3; OR2T34; Olfactory receptor 2T34; Olfactory receptor OR1-63

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2T3/34.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A