

SH3RF3 Antibody



PACO04114

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

SH3RF3

Uniprot

Q8TEJ3

Synonyms:

SH3RF3; POSH2; SH3MD4; SH3 domain-containing RING finger protein 3; Plenty of SH3s 2; SH3 multiple domains protein 4

Immunogen:

Synthesized peptide derived from the C-terminal region of human SH3RF3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A