

# GPR149 Antibody



PACO04115

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

WB:1:500-1:2000, IHC:1:100-1:300,  
IF:1:200-1:1000

**Protein Background:****Gene ID:**

GPR149

**Uniprot**

Q86SP6

**Synonyms:**

GPR149; PGR10; Probable G-protein coupled receptor 149; G-protein coupled receptor PGR10

**Immunogen:**

Synthesized peptide derived from the Internal region of human GPR149.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A