

Phospho-IFNGR1 (Y457) Antibody



PACO04119

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

IFNGR1

Uniprot

P15260

Synonyms:

IFNGR1; Interferon gamma receptor 1; IFN-gamma receptor 1; IFN-gamma-R1;
CDw119; CD antigen CD119

Immunogen:

Synthesized peptide derived from human IFN-gammaRalpha around the phosphorylation site of Y457.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A