

Phospho-IL9R (S519) Antibody



PACO04156

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:

Gene ID:

IL9R

Uniprot

Q01113

Synonyms:

IL9R; Interleukin-9 receptor; IL-9 receptor; IL-9R; CD antigen CD129

Immunogen:

Synthesized peptide derived from human IL-9R around the phosphorylation site of S519.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A