

AQP3 Antibody



PACO04163

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

AQP3

Uniprot

Q92482

Synonyms:

AQP3; Aquaporin-3; AQP-3; Aquaglyceroporin-3

Immunogen:

Synthesized peptide derived from the Internal region of human AQP3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A