

KCNQ2/KCNQ3/KCNQ4/KCNQ5 Antibody



PACO04167

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:**Gene ID:**

KCNQ2/KCNQ3/KCNQ4/KCNQ5

Uniprot

O43526/O43525/P56696/Q9NR82

Synonyms:

KCNQ2; Potassium voltage-gated channel subfamily KQT member 2; KQT-like 2; Neuroblastoma-specific potassium channel subunit alpha KvLQT2; Voltage-gated potassium channel subunit Kv7.2; KCNQ3; Potassium voltage-gated channel subfamily KQT me

Immunogen:

Synthesized peptide derived from human KCNQ2/3/4/5 around the non-phosphorylation site of T217/246/223/251.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A