

OR56A3 Antibody



PACO04189

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR56A3

Uniprot

Q8NH54

Synonyms:

OR56A3; OR56A3P; OR56A6; Olfactory receptor 56A3; Olfactory receptor 56A6

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 56A3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A