

OR5D16 Antibody



PACO04191

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR5D16

Uniprot

Q8NGK9

Synonyms:

OR5D16; Olfactory receptor 5D16; Olfactory receptor OR11-154

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5D16.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A