

OR8K1 Antibody



PACO04195

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

OR8K1

Uniprot

Q8NGG5

Synonyms:

OR8K1; Olfactory receptor 8K1; Olfactory receptor OR11-182

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 8K1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A