

OR5M1/OR5M10 Antibody



PACO04198

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR5M1/OR5M10

Uniprot

Q8NGP8/Q6IEU7

Synonyms:

OR5M1; Olfactory receptor 5M1; OST050; Olfactory receptor OR11-208; OR5M10; Olfactory receptor 5M10; Olfactory receptor OR11-207

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5M1/5M10.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A