

GPR152 Antibody



PACO04200

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

IF:1:200-1:1000

Protein Background:

Gene ID:

GPR152

Uniprot

Q8TDT2

Synonyms:

GPR152; PGR5; Probable G-protein coupled receptor 152; G-protein coupled receptor PGR5

Immunogen:

Synthesized peptide derived from the Internal region of human GPR152.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A