

OR5AU1 Antibody



PACO04204

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR5AU1

Uniprot

Q8NGC0/Q6IEU2

Synonyms:

Olfactory receptor 5AU1

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 5AU1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A