

# OR2AK2 Antibody



PACO04212

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR2AK2

**Uniprot**

Q8NG84

**Synonyms:**

OR2AK2; OR2AK1P; Olfactory receptor 2AK2; Olfactory receptor 2AK1; Olfactory receptor OR1-47

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2AK2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A