

LEP Antibody



PACO04230

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

LEP

Uniprot

P41159

Synonyms:

LEP; OB; OBS; Leptin; Obese protein; Obesity factor

Immunogen:

Synthesized peptide derived from the N-terminal region of human Ob.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A