

Phospho-LIMK2 (S283) Antibody



PACO04238

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

LIMK2

Uniprot

P53671

Synonyms:

LIMK2; LIM domain kinase 2; LIMK-2

Immunogen:

Synthesized peptide derived from human LIMK-2 around the phosphorylation site of S283.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A