

RHOH Antibody



PACO04239

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

RHOH

Uniprot

Q15669

Synonyms:

RHOH; ARHH; TTF; Rho-related GTP-binding protein RhoH; GTP-binding protein TTF;
Translocation three four protein

Immunogen:

Synthesized peptide derived from the C-terminal region of human Rho H.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A