

OR51H1P Antibody



PACO04243

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR51H1

Uniprot

Q8NH63

Synonyms:

OR51H1P; OR51H1; Putative olfactory receptor 51H1; Olfactory receptor OR11-25

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 51H1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A