

SPINK6 Antibody



PACO04253

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

SPINK6

Uniprot

Q6UWN8

Synonyms:

SPINK6; Serine protease inhibitor Kazal-type 6; Kallikrein inhibitor

Immunogen:

Synthesized peptide derived from the Internal region of human SPINK6.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A