

LTK Antibody



PACO04255

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

LTK

Uniprot

P29376

Synonyms:

LTK; TYK1; Leukocyte tyrosine kinase receptor; Protein tyrosine kinase 1

Immunogen:

Synthesized peptide derived from the Internal region of human LTK.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A