

MAP2 Antibody



PACO04278

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MAP2

Uniprot

P11137

Synonyms:

MAP2; Microtubule-associated protein 2; MAP-2

Immunogen:

Synthesized peptide derived from the N-terminal region of human MAP-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A