

**Phospho-MARK1/MARK2/MARK3/MARK4 (T215)
Antibody
PACO04294**



Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

IHC:1:100-1:300

Protein Background:

Gene ID:

MARK1/MARK2/MARK3/MARK4

Uniprot

Q9P0L2/Q7KZI7/P27448/Q96L34

Synonyms:

MARK1; KIAA1477; MARK; Serine/threonine-protein kinase MARK1; MAP/microtubule affinity-regulating kinase 1; PAR1 homolog c; Par-1c; Par1c; MARK2; EMK1; Serine/threonine-protein kinase MARK2; ELKL motif kinase 1; EMK-1; MAP/microtubule affin

Immunogen:

Synthesized peptide derived from human MARK1/2/3/4 around the phosphorylation site of T215.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A