

Phospho-MDM2 (S166) Antibody



PACO04308

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300,
IF:1:200-1:1000

Protein Background:**Gene ID:**

MDM2

Uniprot

Q00987

Synonyms:

MDM2; E3 ubiquitin-protein ligase Mdm2; Double minute 2 protein; Hdm2;
Oncoprotein Mdm2; p53-binding protein Mdm2

Immunogen:

Synthesized peptide derived from human MDM2 around the phosphorylation site of S166.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A