

MDM2 Antibody



PACO04310

Product Information

Size:

50ug

Gene ID:

MDM2

Reactivity:

Human, Mouse

Uniprot

Q00987

Source:

Rabbit

Synonyms:

MDM2; E3 ubiquitin-protein ligase Mdm2; Double minute 2 protein; Hdm2; Oncoprotein Mdm2; p53-binding protein Mdm2

Isotype:

IgG

Immunogen:

Synthesized peptide derived from the C-terminal region of human Mdm2.

Form:

Liquid

Immunogen Species:

Homo sapiens (Human)

Conjugate:

Non-conjugated

Storage Buffer:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Type:

Polyclonal Antibody

Storage:

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Applications:

ELISA, WB, IHC, IF

Purification Method:

The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Recommended dilutions:

WB:1:500-1:2000,

IHC:1:100-1:300,

IF:1:200-1:1000

ELISA:1:40000