

Phospho-MEF2C (S396) Antibody



PACO04316

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

MEF2C

Uniprot

Q06413

Synonyms:

MEF2C; Myocyte-specific enhancer factor 2C

Immunogen:

Synthesized peptide derived from human MEF-2C around the phosphorylation site of S396.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A