

OR10C1 Antibody



PACO04357

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

WB:1:500-1:2000

Protein Background:

Gene ID:

OR10C1

Uniprot

Q96KK4/Q5SUN7

Synonyms:

Olfactory receptor 10C1

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 10C1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A